

Troubleshooting System Blockage

Use this procedure if you have isolated a chromatography/spectrum problem to a restriction or clog in your ACQUITY UPC²® system.

CAUTION: *OBSERVE ALL APPROPRIATE SAFETY AND HANDLING PRECAUTIONS (SEE THE RELEVANT [Overview and Maintenance Guide](#)). IF YOU HAVE ANY DOUBT ABOUT YOUR ABILITY TO PERFORM THIS PROCEDURE, CONTACT WATERS TECHNICAL SERVICE FOR ASSISTANCE.*

Parts required

Table 1: Parts Required^a

Part number	SAP description
700002747	Pin Plug, 1/16", High Pressure, (Pkg/5)
WAT084078	Cap, FM 05 Hex SS
700002636	Union, .020 ID, V-Detail, (1ea)

- a. To ensure that part numbers are up to date, check the appropriate part locator in the Waters Portal.

Disconnecting fittings to isolate blockage

- From the plot properties window, set Amount of Data Shown to **1 minute**.
- In the BSM plot properties window, deselect all plots except for Primary B Pressure and Accumulator B Pressure.
- Remove the CO₂ inlet (the tube coming from P5 on vent valve) from the 250- μ L mixer and install a V-detail pin plug into the mixer inlet.
- Connect a cap plug (WAT084078) to the CO₂ mixer inlet tube.
- Replace the column with a union.
- Run 100% methanol (MeOH) at 1.0 mL/min, with ABPR set to **Off**, Column Position set to location of union, and Column Temp set to **Off**. Allow the accumulator pressure to equilibrate.
- Disconnect fitting at the ABPR inlet (position 1 on [Figure 1](#)).
 - If pressure on the B-side accumulator trace is less than (<) 300 psi, you do not have a blockage before the ABPR.

NOTE: *Do not observe system pressure for this step.*

- If pressure on the B-side accumulator trace is above (>) 300 psi, disconnect fittings, proceeding upstream until the pressure drops below 300 psi. Proceed in the order indicated in [Figure 1](#) and [Table 2](#).

NOTE: *To relieve pressure, slowly loosen the fitting. Use a tissue to prevent solvent spray and to collect any spilled solvent.*

- Once the pressure drops below 300 psi, replace the component downstream from the fitting you just loosened.




Figure 1 - Disconnecting fittings to isolate a blockage

Table 2: Fitting locations

Location number	Description
1	ABPR inlet
2	Transducer outlet
3	Transducer inlet
4	Vent valve tee outlet
5	Vent valve tee inlet
6	PDA outlet
7	PDA inlet
8	Column manager outlet valve center port
9	Column manager outlet valve port 3
10	Column outlet
11	Column inlet or guard column
12	Column manager inlet valve port 3
13	Column manager inlet valve center port
14	Convergence manager auxiliary valve port 6
15	Convergence manager auxiliary valve port 1
16	Sample manager injection valve port 6
17	Sample manager injection valve port 5
18	Convergence manager auxiliary valve port 4
19	Convergence manager auxiliary valve port 15
20	BSM 250- μ L mixer inlet
21	BSM 250- μ L mixer outlet