


# WATERS ION ANALYSIS METHOD

## Anion Analysis Using IC-Pak A HC Column Borate/Gluconate Eluent

Method #  
A-102


### INSTRUMENTATION



### CONDITIONS

**Eluent:** Borate/Gluconate  
**Pump:** 590 Solvent Delivery Module  
**Injector:** 710B WISP  
**Column:** IC-Pak A HC  
**Data:** 840 Data System  
**Flow rate:** 2.0 ml/min  
**Injection:** 100 µl of Standard Anion Mix

**Detection:** 430 Conductivity  
**Range:** 500 µS  
**Temperature:** On  
**Polarity:** +  
**Background:** 274 µS



## ELUENT PREPARATION

### I. Borate/Gluconate Concentrate

To a one liter volumetric flask add:

16 g Sodium Gluconate  
18 g Boric Acid  
25 g Sodium Tetraborate Decahydrate

Add approximately 500 ml of Milli-Q water and mix thoroughly until dissolved, then add:

250 ml Glycerin

Fill the flask to the mark with Milli-Q water and mix thoroughly. Concentrate may be stored refrigerated for up to six months before replacement.

### II. Borate/Gluconate Eluent (pH 8.5)

Place approximately 500 ml of Milli-Q water into a one liter volumetric flask and add:

20 ml Borate/Gluconate Concentrate  
20 ml n-Butanol  
120 ml Acetonitrile

Fill the flask to the mark with Milli-Q water and mix thoroughly. Filter through a 0.22  $\mu\text{m}$  Durapore membrane (GVWP).

---

## STANDARD PREPARATION

### I. Concentrated Standard

1000 ppm $\text{F}^-$	0.221 g NaF (ACS) +/- 0.001 g
2000 ppm $\text{Cl}^-$	0.329 g NaCl (ACS) +/- 0.001 g
4000 ppm $\text{NO}_2^-$	0.600 g $\text{NaNO}_2$ (ACS) +/- 0.001 g
4000 ppm $\text{Br}^-$	0.596 g KBr (ACS) +/- 0.001 g
4000 ppm $\text{NO}_3^-$	0.548 g $\text{NaNO}_3$ (ACS) +/- 0.001 g
6000 ppm $\text{HPO}_4^{2-}$	0.854 g $\text{KH}_2\text{PO}_4$ (ACS) +/- 0.001 g
4000 ppm $\text{SO}_4^{2-}$	0.592 g $\text{Na}_2\text{SO}_4$ (ACS) +/- 0.001 g

Prepare separate standard concentrates by diluting each of the above to 100 ml with Milli-Q water

### II. Working Standard

Pipet 100  $\mu\text{l}$  of each standard concentrate into a 100-ml volumetric flask and dilute with Milli-Q water to result in the working standard (see front page). Prepare fresh working standard weekly.

**DETECTION LIMITS:**

3 x Baseline Noise (500 µl)

<u>Analyte</u>	<u>PPB</u>
F <sup>-</sup>	20
Cl <sup>-</sup>	10
NO <sub>2</sub> <sup>-</sup>	20
Br <sup>-</sup>	30
NO <sub>3</sub> <sup>-</sup>	20
HPO <sub>4</sub> <sup>2-</sup>	75
SO <sub>4</sub> <sup>2-</sup>	50

**COMMENTS:**

1. For samples with low levels of analytes, injection volumes of 500 µl may be used.
2. For specific sample preparation methods, see Section 1.
3. When analyzing for anions in drinking water and other environmental samples, the presence of high levels of carbonate, interfering metals (such as magnesium and calcium), and excess base are often encountered. In an untreated sample, these may cause disturbances in the chromatographic baseline. The presence of high levels of carbonate will result in a large carbonate peak which may interfere with the quantitation of fluoride and chloride. In extreme cases carbonate may overload the column and cause loss of early eluting analytes. Alkaline earth metals, such as magnesium and calcium, may cause severe interferences in the early part of the chromatogram.

When using a borate/gluconate eluent, high pH samples may result in a pH system peak in the middle portion of the chromatogram which can interfere with nitrate and bromide analysis. In many of these samples, the analyte levels are too low for dilution of the sample to be practical. Dilution would cause analyte levels to drop below the detection limits for direct injection.

The best method for analyzing samples of this type is to pretreat the sample with a MilliTrap™ H<sup>+</sup> Membrane Cartridge. The MilliTrap™ H<sup>+</sup> cartridge is a hand held, multiple use, disposable device that removes interferences from many matrices such as drinking, surface or waste water samples. The diagram on pages 7 & 8 of Section 1 illustrates and explains how the MilliTrap™ H<sup>+</sup> cartridge functions to remove matrix interferences.

When using the MilliTrap™ H+ Membrane Cartridge, the eluent described below is recommended. Using the MilliTrap™ cartridge with this eluent will give superior results for both high pH and low pH samples.

---

### Concentrate

To a 1-liter volumetric flask add:

500 ml Milli-Q® water

7.2 g Lithium Hydroxide Monohydrate

Mix until dissolved, then add:

25.5 g Boric Acid

Mix until dissolved, then add:

13.2 ml Gluconic Acid (50% wt solution)

94 ml Glycerin (95%)

Finally, add Milli-Q® water to the mark and mix thoroughly.

### Eluent

To a 1-liter volumetric flask add:

20 ml of above concentrate

120 ml Acetonitrile

Fill to the mark with Milli-Q® water, mix thoroughly, and filter through a 0.22 µ Millipore membrane (GVWP). Background conductivity = 220 µS.

---